Xiamen Boson Biotech Co.,Ltd.

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www.bosonbio.com

TABLE OF CONTENTS

1	GENERAL APPLICATION INFORMATION	4
2	PRODUCT DESCRIPTION	5
2.1	General Description	5
2.2	Intended Use & Operation.....	5
2.2.1	Diagnostic Principle	5
2.2.2	Intended Use	5
3	PHOTOGRAPHS OF COMPONENTS.....	6
3.1	Comprehensive Description of the Product.....	6
3.2	Photographs of Components	7
4	VALIDATION STUDIES.....	10
4.1	Analytical Performance Studies of the Test in Professional Use.....	10
4.2	Stability Studies	10
4.3	Clinical Evidence.....	10
4.4	Analytical Performance Studies of the Test in Self-Use	10
4.5	User-Friendliness Study for Self-Test	10
5	LABELS & INSTRUCTIONS FOR USE.....	12
5.1	Labels	12
5.1.1	For 1 Test/Box.....	12
5.1.2	For 5 Tests/Box	13
5.2	Instructions for Use (IFU).....	15
5.2.1	Instructions for Use in English	15
5.2.2	Instructions for Use in Dutch	20
6	SPECIFICATION OF MATERIALS.....	25
6.1	Compositions of the Test Kit	25
7	CERTIFICATIONS	26
7.1	ISO13485 Certificate.....	26
7.2	CE Certificate of Professional Use Test.....	28
7.3	Declaration of Conformity of Professional Test	30
7.4	CE Certificate of Other Components (N/A)	31
7.5	EU Authorization	31

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Tel : 5 1 2e E-mail: 5 1 2e

www.bosonbio.com

7.6	BfArM Authorization for Self-Test.....	32
7.7	Confirmation of CE Certificate Application for Use as Self-Test.....	32
8	ESSENTIAL REQUIREMENTS CHECKLIST.....	34
9	RISK MANAGEMENT.....	50
10	LIST OF APPENDICES.....	52

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www.bosonbio.com

1 GENERAL APPLICATION INFORMATION

This application is submitted on behalf of Xiamen Boson Biotech Co., Ltd., and is intended as an exemption for placing the Rapid SARS-CoV-2 Antigen Test Card on the Dutch market as a self-test. Current regulatory approvals and processes are listed below:

1. The Rapid SARS-CoV-2 Antigen Test Card manufactured by Xiamen Boson Biotech Co., Ltd. has received the CE mark for professional use as an antigen rapid test (See Section 7.2 and 7.3 for the CE certificate and Declaration of Conformity).
2. We have already started the conformity assessment procedure to obtain a CE certificate for the use of this antigen rapid test as a self-test via TÜV SÜD. See Section 7.7 for the confirmation email for the CE submission).
3. The professional test is included in the most recent version of the op 17 February 2021 adopted list of the Health Security Committee "A common list of COVID-19 rapid antigen tests". See Section 7.5 for the screenshot of the document with the listed product, and check the following website for the official document:

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4. The home test is special approved by BfArM in Germany as a test for personal use by laypeople. See Section 7.6 for the screen shot from the official website of BfArM and check the following website for the original contents:

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2 PRODUCT DESCRIPTION

2.1 General Description

Product / Device name:	Rapid SARS-CoV-2 Antigen Test Card
Catalog Number:	1N40C5-2 (1 Test/Box) 1N40C5-4 (5 Tests/Box)
Mode of Action	RDT (Rapid diagnostic test)

2.2 Intended Use & Operation

2.2.1 Diagnostic Principle

Rapid SARS-CoV-2 Antigen Test Card is an immunochromatographic lateral flow device that employs the principle of double antibody sandwich method. Colloidal gold conjugated anti-SARS-CoV-2 antibodies are dry-immobilized on the test device. When the specimen is added, it migrates by capillary diffusion through the strip to re-hydrate the gold conjugate complexes. If present at or above the limit of detection, SARS-CoV-2 viral antigens will react with the gold conjugate complexes to form particles, which will continue to migrate along the strip until the Test Zone (T) where they are captured by the immobilized anti-SARS-CoV-2 antibodies to form a visible red line. If there are no SARS-CoV-2 viral antigens in the specimen, no red line will appear in the Test Zone (T). The gold conjugate complexes will continue to migrate alone until being captured by immobilized antibody in the Control Zone (C) to form a red line, which indicates the validity of the test.

2.2.2 Intended Use

Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card cannot be used as the basis to diagnose or exclude SARS-CoV-2 infection.

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 2e E-mail: 5 1 2e
www.bosonbio.com

3 PHOTOGRAPHS OF COMPONENTS

3.1 Comprehensive Description of the Product

The Rapid SARS-CoV-2 Antigen Test Card is a cassette device for detecting the presence or absence of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. The plastic housing of device is designed to affix the test strip inside the device. There is no functional design that may affect the test strip's performance in the test card.

Each test strip contains an absorbent pad, coated membrane, gold conjugate pad, and sample pad with adhesive sticker. The test strip is used to assemble the test card device.

Components	For 1 Test/Box	For 5 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5
Sterilized swab	1	5
Extraction tube	1	5
Sample extraction buffer	1	5
Instructions for use	1	1
Tube stand	1 (In Box)	1

For 1 Test/Box: All components of the test kit are the same to the original professional use kit.

For 5 Tests/Box: All components of the test kit are the same to the original professional use kit, except for the different sample extraction buffer specifications. For the professional use kit, the sample extraction buffer is 1 bottle for 2mL. For the self-test kit, the sample extraction buffer is 5 single packs.

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www.bosonbio.com

3.2 Photographs of Components

Picture of all components for 1 test/kit

Picture of all components for 5 tests/kit

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
Tel : 5 1 2e E-mail: 5 1 2e

www.bosonbio.com


Picture of swab**Picture of sample extraction buffer****Picture of extraction tube**

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 2e E-mail: 5 1 2e
www.bosonbio.com

Picture of test card

4 VALIDATION STUDIES

4.1 Analytical Performance Studies of the Test in Professional Use

No.	Study	Reference
1	Analytical Sensitivity	See Appendix 3-1 - Analytical Sensitivity Studies See Appendix 3-2 - Analytical Sensitivity Studies- Nasal Sample Validation Report
2	Analytical Specificity- Cross Reactivity	See Appendix 4 – Analytical Specificity – Cross Reactivity Studies
3	Analytical Specificity- Interference	See Appendix 5 – Analytical Specificity – Interference Studies
4	Precision Studies	See Appendix 6 – Precision Studies
5	High Dose Hook Effect	See Appendix 7 – Hook Effect Studies

4.2 Stability Studies

No.	Study	Reference
1	Storage Stability	See Appendix 8 – Storage Stability Studies
2	Freeze-Thaw Stability	See Appendix 9 – Freeze-Thaw Stability Studies
3	Transport Stability	See Appendix 10 – Transport Stability Studies
4	Heat Stability	See Appendix 11 – Heat Stability Studies
5	Open-Pack Stability	See Appendix 12 – Open Pack Studies

4.3 Clinical Evidence

See Appendix 13 – Clinical Evaluation Report.

4.4 Analytical Performance Studies of the Test in Self-Use

Because the components for the professional and self-use tests are identical, the analytical performance studies of the test in self-use are the same as for the professional test.

4.5 User-Friendliness Study for Self-Test

See Appendix 14 – Private-Use Usability Study.

The private use usability study was completed in Austria. Considering that the majority of

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Tel : [REDACTED] E-mail: [REDACTED]

www.bosonbio.com

population in both the Netherlands and Austria are of the Caucasian origin, it is safe to extrapolate the test results in the report to the Dutch situation.

In the future, we also plan to carry out more layperson studies and validate the self-tests in the Netherlands. We will submit follow-up reports once they are available.

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www.bosonbio.com

5 LABELS & INSTRUCTIONS FOR USE

The requirements of Annex I, B, 8 ff in Directive 98/79/EC are fulfilled, as shown on the following pages. An easy-to-understand and easy-to-read instruction leaflet for safe and correct use is included with each kit.

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

For private use/home use/self-testing

REF	1N40C5-2	For 1 Test/Box
REF	1N40C5-4	For 5 Tests/Box

Please follow the instruction leaflet carefully.

INTENDED USE

Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

MATERIALS PROVIDED

Components	For 1 Test/Box	For 5 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5
Sterilized swab	1	5
Extraction tube	1	5
Sample extraction buffer	1	5
Instructions for use (this leaflet)	1	1
Tube stand	1 (In box)	1

IMPORTANT INFORMATION BEFORE THE EXECUTION

1. Read this instruction guide carefully.
2. Do not use the product beyond the expiration date.
3. Do not use the product if the pouch is damaged or the seal is broken.
4. Store the test device at 4 to 30°C in the original sealed pouch. Do Not Freeze.
5. The product should be used at room temperature (15°C to 30°C). If the product has been stored in a cool area (less than 15°C), leave it at normal room temperature for 30 minutes before using.
6. Handle all specimens as potentially infectious.
7. Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
8. Use the swabs included in the test kit to ensure optimal performance of the test.
9. Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.
10. Blow the nose several times before collecting specimen.
11. The specimens should be tested as soon as possible after collection.
12. Apply the drops of test specimen only to the specimen well (S).
13. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
14. Children under 14 years of age should be assisted by an adult.

LIMITATIONS

1. The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.
2. Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results. Improper collection, storage,

or even freezing and thawing of the specimen can lead to inaccurate test results.

3. If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.
5. A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.
6. A positive result does not exclude coinfection with other pathogens.
7. The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may not correlate with other diagnostic methods performed on the same specimen.
8. Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.
9. Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.
10. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
11. The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.
12. The kit was validated with the assorted swabs. Use of alternative swabs may result in false negative results.
13. The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.

PREPARATION


- Clear, clean and dry a flat surface.
- Check the test kit contents. Make sure that nothing is damaged or broken.
- Timer at hand.
- Blow your nose several times before collecting specimen.
- Wash hands.

DISPOSAL

The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

PROCEDURE

This test is suitable for people of all ages. The recommended operators are aging from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child feels any pain.

- 1  Rotate the lid of sample extraction buffer bottle.
Caution: Open it away from your face and be careful not to spill any of the liquid.
- 2  Squeeze all extraction buffer out of the bottle into the extraction tube.
Caution: Avoid touching the bottle against the tube.
- 3  Find the swab in the sealed wrapper in front of you. Identify the soft, fabric tip of the swab.
- 4  Peel open the swab packaging and gently take out the swab.
Caution: Never touch the soft, fabric tip of the swab with your hands.


Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
 Xiamen, Fujian 361021, P.R. China
 Tel : 5 1 2e E-mail: 5 1 2e
 www.bosonbio.com


5 Carefully insert swab into one nostril. The swab tip should be inserted no less than 2.5 cm (1 inch) from the edge of the nostril. Roll swab 3-4 times along the mucosa inside the nostril. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril. Withdraw swab from the nasal cavity.


Caution: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.


6 Place swab into extraction tube. Roll swab three to five (3-5) times. **Leave swab in extraction buffer for 1 minute.**

7 Pinch extraction tube with fingers and remove the solution from swab as much as possible.


8 Install the nozzle cap onto the sample extraction tube tightly.

9 Bring the kit components to room temperature before testing. Open the pouch and remove the card. Place the card on a flat and level surface.

Caution: Once opened, the test card must be used immediately.

10 Invert the extraction tube and add 3 drops of test specimen into the specimen well (S) by gently squeezing the extraction tube.

Caution: The formation of air bubbles in the specimen well (S) must be avoided.

11 Read the results at 15-20 minutes.

Caution: Results after 20 minutes may not be accurate. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

INTERPRETATION OF RESULTS

Positive: If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.

Caution: No matter how faint the colored band is in the Test Zone (T), the result

should be considered as positive.

Negative:

If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.

Invalid:

If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

QUALITY CONTROL

The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagents are reactive.

FREQUENTLY ASKED QUESTIONS (FAQ)
1. How does the detection work?

The N protein of the SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears. Therefore, if the sample does not contain any viral proteins or antigens, there will be no red test line (T).

2. When should/can I test myself?

You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.

3. What can affect my test result? What should I pay attention to?

Be sure to blow your nose multiple times before collecting the specimen.
 Be sure to visibly collect sample material (nasal secretions).
 Perform the test immediately after taking the sample.
 Follow the instructions for use carefully.
 Apply the drops of extraction solution only to the sample well (S)
Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.

4. The test strip is clearly discolored or smudged? What is the reason for this?

Please note that the test card should not be used with more than 3 drops of sample, as the liquid absorption of the test strip is naturally limited. If the control line does not appear or the test strip is badly smudged or discolored, making it unreadable, please repeat the test according to the instructions.

5. I have taken the test, but I don't see a control line (C). What should I do?

Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.

6. I am unsure about reading the result. What should I do?

For the result to be positive, 2 straight horizontal lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest health facility according to the recommendations of your local authorities.

7. My result is positive. What should I do?

If your result is positive and the test kit thus clearly indicates the control line as well as the test line, please contact your local GGD (municipal health service) for further instructions. Your test result may be double-checked and the GGD will explain the appropriate next steps.

8. My result is negative. What should I do?

If the test kit only clearly shows the control line, this may mean that you are negative or that the viral load is too low to be detected. If you experience symptoms (headache, fever, migraine, loss of sense of smell or taste, etc.), please consult your primary care physician, or the nearest health care facility as recommended by your local authorities.
 If the test results are consistently negative for 5 days after symptom onset, you can stop self-quarantining.
 If you are still unsure about the test results, make an appointment via www.government.nl/coronavirus-test, call 0800 1202, or contact your local GGD for further instructions.

9. How can I dispose of the product?


The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

BOSON

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 2e E-mail: 5 1 2e
www.bosonbio.com

EXPLANATION FOR SYMBOLS

	In Vitro Diagnostics Use		See Instructions for Use		Expiry Date
	Tests per Kit		Keep dry		Batch Number
	Authorized Representative		Keep away from sunlight		Manufacturer
	Do not reuse		Do not use if package is damaged		Store between 4- 30 °C
	CE Mark		Catalogue Number		Warning, please refer to the instruction

Manufacturer:  Xiamen Boson Biotech Co., Ltd.
90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian, 361021, P.R.China.

Authorized Representative: Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands.

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6 SPECIFICATION OF MATERIALS

6.1 Compositions of the Test Kit

Component	Amount	Specification	Function	Supplier
Conjugate Pad	0.2 µg/test unit	Au-Mouse anti-SARS-CoV-2 antibodies Conjugate	Conjugating antibody	Hangzhou Clongene Biotech Co., Ltd.
Test membrane, test line	0.6 µg/test unit	Mouse anti-SARS-CoV-2 antibodies	Coating antibody	Hangzhou Clongene Biotech Co., Ltd.
Test membrane, control line	0.6 µg/test unit	Goat anti mouse IgG antibody	Quality control	Boson
Nitrocellulose Membrane	1.83 cm × 0.4 cm	Nitrocellulose	Coating control line and test line	Shenzhen Bai Sui Kang Industrial Co., Ltd.
Glass Fiber Paper	1.7 cm × 0.4 cm	Glass Fiber	Sample addition	Hangzhou Huarun Paper Co., Ltd.
Absorbent Paper	2.27 cm × 0.4 cm	Wood Fiber	Keeping capillary attraction	Hangzhou Special Paper Industry Co., Ltd. (Newstar)
Swab	1 pc or 5 pcs	Flocking and PP	Collecting nasal samples	Goodwood Medical Care Ltd.
Sample Extraction Buffer	0.27 mL/Bottle, 1 bottle or 5 bottles	Na ₂ HPO ₄ ·12H ₂ O, NaH ₂ PO ₄ ·2H ₂ O, NaCl, TERGITOL 15-S-9 and Proclin300	Lysis virus	Boson

All the specifications, functions and suppliers of the components for the self-test are the same as for the professional use test.

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Certificate

No. Q5 061317 0005 Rev. 00

Holder of Certificate: **Xiamen Boson Biotech Co., Ltd.**
90-94 Tianfeng Road
Jimei North Industrial Park
361021 Xiamen, Fujian
PEOPLE'S REPUBLIC OF CHINA

Facility(ies): Xiamen Boson Biotech Co., Ltd.
90-94 Tianfeng Road, Jimei North Industrial Park, 361021
Xiamen, Fujian, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

tuv-sud.com/pe-cert

Scope of Certificate: Design and Development, Production and Distribution of
In Vitro Diagnostic kits for detection of Infection Diseases,
Tumour Markers, Drug Abuse, Hormones, Cardiac Markers
and Related Biomaterial

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH1807513

Valid from: 2018-10-31

Valid until: 2021-10-30

5.1.2e

Date, 2018-09-27

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,

Xiamen, Fujian 361021, P.R. China

Tel : 5 1 2e E-mail: 5 1 2e

www.bosonbio.com

Smart Certificate Explorer

Your scan returns the following certificate:

No.
Q5 061317 0005 Rev. 00

QM Certificates
QM System ISO 13485

Holder of Certificate:
Xiamen Boson Biotech Co., Ltd.
90-94 Tianfeng Road
Jimei North Industrial Park
361021 Xiamen, Fujian
PEOPLE'S REPUBLIC OF CHINA

-

Scope of Certificate:
Design and Development, Production and Distribution of In Vitro Diagnostic kits for detection of Infection Diseases, Tumour Markers, Drug Abuse, Hormones, Cardiac Markers and Related Biomaterial

-

Date.
2018-09-27

Valid

-

-

-

-

-

Note: Valid certificates are reported only. Not all information is available in English. Binding information and information about withdrawn certificates can be obtained from our certification department
512e@tuv-sued.de
[→ Link to Certificate Explorer](#)

[Imprint](#) [Privacy Statement](#)

7.2 CE Certificate of Professional Use Test

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Farmatec

[REDACTED] 5.1.2e

7.3 Declaration of Conformity of Professional Test

 Declaration of Conformity	
Manufacturer	Xiamen Boson Biotech CO., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P. R. China.
European Representative	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Device/s	Rapid SARS-CoV-2 Antigen Test Card
Classification	Others
Confirmative Assessment Route	98/79/EC IVDD Annex III
<p>We, Xiamen Boson Biotech Co., Ltd. declare that the above mentioned devices conforms to the relevant provisions of the EC Council Directive 98/79 and is in accordance with the Annex III, ISO 13485:2016 Quality Management System, as implemented by the European Union's Medical Devices Regulations and the Federal and Local Authorities.</p>	
Place, Date of Issue	Xiamen, 2020-08-04
Signature	[redacted] 5.1.2e (Signed By Boson Representative) Name: [redacted] 5.1.2e Title: [redacted] 5.1.2e
	

7.4 CE Certificate of Other Components (N/A)

Not applicable. The components of the self-test is the same as for the professional test.

7.5 EU Authorization

The professional test is included in the most recent version of the op 17 February 2021 adopted list of the Health Security Committee "A common list of COVID-19 rapid antigen tests". Below is the screenshot taken from page 4 of the original document. Refer to the following link for the full official document:

https://ec.europa.eu/health/sites/health/files/preparedness_response/docs/covid-19_rat_common-list_en.pdf

information in this platform and ensuring that the common list as agreed by Member States will be publicly available online.

The common list of rapid antigen tests will be regularly reviewed by Member States in the context of Health Security Committee meetings, and, if necessary, be updated in line with new results from independent validation studies becoming available and new tests entering the markets. Future updates to the list should also take into account how mutations of the SARS-CoV-2 virus may affect the efficacy of any particular rapid antigen tests, allowing for the removal of tests no longer deemed effective. The effect of mutations of the SARS-CoV-2 virus on the efficacy of NAAT, in particular RT-PCR assays, will also be kept under review.

Future updates to the common list of rapid antigen tests will be published as an update to the JRC database on COVID-19 In Vitro Diagnostic Devices and Test Methods.

III. Rapid antigen tests of which the test results are mutually recognised

As stipulated in point 15 of the Council Recommendation of 21 January 2021, Member States will agree on a selection of rapid antigen tests of which they will mutually recognise the test results for public health measures, based on the information included in the common list (see Annex I).

The Health Security Committee agrees that, for rapid antigen test results to be mutually recognised, at least three Member States should be using a rapid antigen tests in practice. Based on this criterion, Member States agree that the results of the following rapid antigen tests will be mutually recognised for public health measures:

- Abbott Rapid Diagnostics, Panbio™ COVID-19 Ag Rapid Test
- AMEDA Labordiagnostik GmbH, AMP Rapid Test SARS-CoV-2 Ag
- Becton Dickinson, BD Veritor System for Rapid Detection of SARS-CoV-2
- Beijing Lepu Medical Technology, SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold immunochromatography)
- BIOSYNEX SWISS SA, BIOSYNEX COVID-19 Ag BSS
- CerTest Biotech S.L., CerTest SARS-CoV-2 CARD TEST
- Hangzhou Clongene Biotech, Clungene COVID-19 Antigen Rapid Test Kit
- Healgen Scientific Limited, Coronavirus Ag Rapid Test Cassette (Swab)
- LumiraDX UK LTD, LumiraDx SARS-CoV-2 Ag Test
- nal von minden GmbH, NADAL COVID -19 Ag Test
- Quidel Corporation, Sofia 2 SARS Antigen FIA
- SD BIOSENSOR, Inc., STANDARD F COVID-19 Ag FIA
- SD BIOSENSOR, Inc., STANDARD Q COVID-19 Ag Test
- Siemens Healthineers, CLINITEST Rapid COVID-19 Antigen Test
- **Xiamen Boson Biotech Co, Rapid SARS-CoV-2 Antigen Test card**
- Zhejiang Orient Gene Biotech Co.,Ltd, Coronavirus Ag Rapid Test Cassette (Swab)

7.6 BfArM Authorization for Self-Test

The home test is special approved by BfArM in Germany as a test for personal use by laypeople. Below is the screen shot from the official website of BfArM. Refer to the following link for the original contents:

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7.7 Confirmation of CE Certificate Application for Use as Self-Test

See below for a confirmation email sent by TÜV SÜD on confirmation of CE certificate application for use as self-test.

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : [REDACTED] E-mail: [REDACTED]
www.bosonbio.com

acceptance of the application for COVID-19 self testing

Weng, Chenchuan <[REDACTED]@tuvsud.com>

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Dear Mr. Haolong Shen,

Greeting from TÜV SÜD, notified body of IVDD 98/79/EC.

This is to confirm that the manufacturer

Xiamen Boson Biotech Co., Ltd., 90-94 Tianfeng Road, Jimei North Industrial Park, 361021 Xiamen, Fujian, People's Republic of China
have initiated the application process for a SARS-CoV-2 Rapid Antigen Test intended for self-testing acc. to IVDD (98/79/EC), Annex III.6. TÜV SÜD has already given the first feedback to Xiamen Boson.

Please note: the acceptance of the application does not guarantee a successful certification of the product.

Kind regards,

[REDACTED]

[REDACTED] In Vitro Diagnostic Medical Device

PS:MHS

Phone: +86 [REDACTED]

Fax: +86 [REDACTED]

Email address: [REDACTED]@tuvsud.com

<http://www.tuv-sud.cn>

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

8 ESSENTIAL REQUIREMENTS CHECKLIST

Device: Rapid SARS-CoV-2 Antigen Test Card

Applicable REF:

1N40C5 Rapid SARS-CoV-2 Antigen Test Card

Manufacturer: Xiamen Boson Biotech Co., Ltd.

Date: 2021.02.06

A. GENERAL REQUIREMENTS

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
1.	The device must be designed and manufactured in such a way that, when used under the conditions and for the purposes intended, they will not compromise, directly or indirectly, the clinical condition or the safety of patients, or the safety and health of users or, where applicable, other persons, or the safety of property. Any risks which may be associated with their use must be acceptable when weighed against the benefits to the patient and be compatible with a high level of protection of health and safety.	Yes	EN ISO 14971:2012 Medical devices - Application of risk management to medical devices.	Product Dossier Section 4.1, Analytical Performance Studies of the Test in Professional Use. Product Dossier Section 4.3, Clinical Evidence. Section 4.4, Analytical Performance Studies of the Test in Self-Use. Product Dossier Section 9, RISK MANAGEMENT	Design and Development Control Procedures BS-QA-008 DM-00-001 Clinical Trial Management Protocol BS-QA-007 Risk Management Control Procedures
2.	The solutions adopted by the manufacturer for the design and construction of the devices must conform to safety principles, taking account of the generally acknowledged state of the art. In selecting the most appropriate solutions, the manufacturer must apply the following principles in the following order: - eliminate or reduce risks as far as possible (inherently safe design and construction)	Yes	EN ISO 14971:2012 Medical devices - Application of risk management to medical devices.	Product Dossier Section 9, RISK MANAGEMENT. Product Dossier Section 5.2, Instructions for Use (IFU). Precautions and Limitations section.	BS-QA-007 Risk Management Control Procedures

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

	- where appropriate take adequate protection measures including alarms if necessary, in relation to risks that cannot be eliminated. - inform users of the residual risks due to shortcomings of the protection measures adopted.				
3.	The devices must be designed and manufactured in such a way that they are suitable for the purposes referred to in Article 1(2)(b), as specified by the manufacturer, taking account of the generally acknowledged state of the art. They must achieve the performances, in particular, where appropriate, in terms of analytical sensitivity, diagnostic sensitivity, analytical specificity, diagnostic specificity, accuracy, repeatability, reproducibility, including control of known relevant interference, and limits of detection, stated by the manufacturer. The traceability of values assigned to calibrators and/or control materials must be assured through available reference measurement procedures and/or available reference materials of a higher order.	Yes	ISO 13485:2016: Medical devices - Quality management systems. EN ISO 13612:2002: Performance evaluation of in vitro diagnostic medical devices. EN ISO14971:2012: Medical devices - Application of risk management to medical devices.	Product Dossier Section 4.1, Analytical Performance Studies of the Test in Professional Use. Product Dossier Section 4.3, Clinical Evidence. Section 4.4, Analytical Performance Studies of the Test in Self-Use. Product Dossier Section 9, RISK MANAGEMENT.	BS-QA-008 Design and Development Control Procedures DM-00-001 Clinical Trial Management Protocol BS-QA-007 Risk Management Control Procedures
4.	The characteristics and performance referred to in sections 1 and 3 must not be adversely affected to such a degree that the health or the safety of the patients or the user and, where applicable, of other persons are compromised during the lifetime of the device as indicated by the manufacturer, when the device is subjected to the stresses which can occur during normal conditions of use. When no lifetime is stated, the same applies for the lifetime reasonably to be expected of a device of that kind, having regard to the intended purpose and the anticipated use of the device.	Yes	EN ISO 13640:2002: Stability testing of in vitro diagnostic medical devices.	Product Dossier Section 4.2, Stability Studies.	Stability Inspection Management Protocol QM-00-032
5.	The devices must be designed, manufactured and packed in such a way that their characteristics and performances during their intended use will not be adversely affected during transport and storage (temperature, humidity, etc.) taking into account the instructions and information provided by the manufacturer.	Yes	ISO 13485:2016 EN ISO 13640:2002: Stability testing of in vitro diagnostic medical devices.	Product Dossier Section 4.2, Stability Studies.	Stability Inspection Management Protocol QM-00-032

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

B. DESIGN AND MANUFACTURING REQUIREMENTS

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
1. 1.1	<u>Chemical and physical properties</u> The devices must be designed and manufactured in such a way as to achieve the characteristics and performances referred to in Section A on the "General Requirements". Particular attention must be paid to the possibility of impairment of analytical performance due to incompatibility between the materials and the specimens (such as biological tissues, cells, body fluids and micro-organisms) intended to be used with the device, taking account of its intended purpose.	Yes	ISO 13485:2016 EN ISO 13612:2002: Performance evaluation of in vitro diagnostic medical devices.	Product Dossier Section 6, SPECIFICATION OF MATERIALS.	BS-QA-008 Design and Development Control Procedures
1.2	The devices must be designed, manufactured and packed in such a way as to reduce as far as possible the risk posed by product leakage, contaminants and residues to the persons involved in the transport, storage and use of the devices, taking account of the intended purpose of the product.	Yes	EN ISO 14971:2012 Medical devices - Application of risk management to medical devices.	Product Dossier Section 9, RISK MANAGEMENT.	BS-QA-007 Risk Management Control Procedures
2. 2.1	<u>Infection and microbial contamination</u> The devices and their manufacturing processes must be designed in such a way as to eliminate or reduce as far as possible the risk of infection to the user or other persons. The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes.	Yes	EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.	Product Dossier Section 5.2, Instructions for Use (IFU). Clinical sample volume and application defined in Procedure section, specimen handling and proper disposal defined in Important Information before the Execution section.	BS-QA-008 Design and Development Control Procedures BS-QA-007 Risk Management Control Procedures

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
2.2	Where a device incorporates biological substances, the risks of infection must be reduced as far as possible by selecting appropriate donors and appropriate substances and by using appropriate validated inactivation, conservation, test and control procedures.	N/A			
2.3	Devices labelled either as 'STERILE' or as having a special microbiological state must be designed, manufactured and packed in an appropriate pack, according to procedures suitable for ensuring that they remain in the appropriate microbiological state indicated on the label when placed on the market, under the storage and transport conditions specified by the manufacturer, until the protective packaging is damaged or opened.	N/A			
2.4	Devices labelled either as 'STERILE' or as having a special microbiological state must have been processed by an appropriate, validated method.	N/A			
2.5	<p>Packaging systems for devices other than those referred to in section 2.3 must keep the product without deterioration at the level of cleanliness indicated by the manufacturer and, if the devices are to be sterilised prior to use, reduce as far as possible the risk of microbial contamination.</p> <p>Steps must be taken to reduce as far as possible microbial contamination during selection and handling of raw materials, manufacture, storage and distribution where the performance of the device can be adversely affected by such contamination.</p>	Yes	ISO 13485:2016	<p>Each test device is packaged in a moisture barrier pouch with 0.5g desiccant. Raw material must meet material specification requirements before release.</p> <p>Production environment for each process specified in the manufacturing process flow chart.</p> <p>Except for box packaging, all production processes are carried in Level 100,000 clean area.</p>	<p>Production Environment Control Procedures BS-QA-010</p> <p>SARS-CoV-2 Antigen Process Protocol PS-11-038</p> <p>Inspection documents and standards for various raw materials QO-</p>

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
2.6	Devices intended to be sterilized, must be manufactured in appropriately controlled (e.g. environmental) conditions.	N/A			
2.7	Packaging systems for non-sterile devices must keep the product without deterioration in the level of cleanliness stipulated and, if the devices are to be sterilized prior to use, minimize the risk of microbial contamination; the packaging system must be suitable taking account of the method of sterilization indicated by the manufacturer.	Yes	ISO 13485:2016 EN ISO 14971:2012 Medical devices - Application of risk management to medical devices.	Production environment for each process specified in the manufacturing process flow chart. Except for box packaging, all production processes are carried in Level 100,000 clean area.	Production Environment Control Procedures BS-QA-010
3 3.1	<u>Manufacturing and environmental properties</u> If the device is intended for use in combination with other devices or equipment, the whole combination, including the connection system must be safe and must not impair the specified performances of the devices. Any restrictions on the use must be indicated on the label or in the instructions for use.	N/A			
3.2	Devices must be designed and manufactured in such a way as to reduce as far as possible the risks linked to their use in conjunction with materials, substances and gases with which they may come into contact during normal conditions of use.	N/A			
3.3	Devices must be designed and manufactured in such a way as to remove or reduce as far as possible: — the risk of injury linked to their physical features (in particular aspects of volume x pressure, dimension and, where appropriate, ergonomic features), — risks linked to reasonably foreseeable external influences, such as magnetic fields, external electrical effects, electrostatic discharge, pressure, humidity, temperature or variations in pressure or acceleration or accidental penetration of substances into the device.	N/A Yes	EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.	Humidity and temperature reference Product Dossier Section 5.2, Instructions for Use (IFU).	Stability Management Protocol QM-00-032

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	Devices must be designed and manufactured in such a way as to provide an adequate level of intrinsic immunity of electromagnetic disturbance to enable them to operate as intended.	N/A	EN ISO 13640: 2002: Stability testing of in vitro diagnostic medical devices.		
3.4	Devices must be designed and manufactured in such a way as to reduce as far as possible the risks of fire or explosion during normal use and in single fault condition. Particular attention must be paid to devices whose intended use includes exposure to or use in association with flammable substances or substances which could cause combustion.	Yes	EN ISO 14971:2012 Medical devices - Application of risk management to medical devices.	Product Dossier Section 9, RISK MANAGEMENT.	BS-QA-007 Risk Management Control Procedures
3.5	Devices must be designed and manufactured in such a way as to facilitate the management of safe waste disposal.	Yes	EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.	Product Dossier Section 5.2, Instructions for Use (IFU). Important Information before the Execution section.	BS-QA-008 Design and Development Control Procedures
3.6	The measuring, monitoring or display scale (including colour change and other visual indicators) must be designed and manufactured in line with ergonomic principles, taking account of the intended purpose of the device.	Yes	EN ISO 13612:2002: Performance evaluation of in vitro diagnostic medical devices.	Device is intended for visual qualitative interpretation.	BS-QA-008 Design and Development Control Procedures BS-QA-007 Risk Management Control Procedures
4 4.1	<u>Devices which are instruments or apparatus with measuring function</u> Devices which are instruments or apparatus having a primary analytical measuring function must be designed and manufactured in such a way as to provide adequate stability and accuracy of measurement within appropriate accuracy limits, taking into account the intended purpose of the device and of available and appropriate reference measurement procedures and materials. The accuracy limits have to be specified by the manufacturer.	N/A			

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
4.2	When values are expressed numerically, they must be given in legal units conforming to the provisions of Council Directive 80/181/EEC of 20 December 1979 on the approximation of the laws of the Member States relating to units of measurement.	N/A			
5	<u>Protection against radiation</u>	N/A			
5.1	Devices shall be designed, manufactured and packaged in such a way that exposure of users and other persons to the emitted radiation is minimised.	N/A			
5.2	When devices are intended to emit potentially hazardous, visible and/or invisible radiation, they must as far as possible be: — designed and manufactured in such a way as to ensure that the characteristics and the quantity of radiation emitted can be controlled and/or adjusted, — fitted with visual displays and/or audible warnings of such emissions.	N/A			
5.3	The operating instructions for devices emitting radiation must give detailed information as to the nature of the emitted radiation, means of protecting the user, and on ways of avoiding misuse and of eliminating the risks inherent in installation.	N/A			
6	<u>Requirements for medical devices connected to or equipped with an energy source</u>	N/A			
6.1	Devices incorporating electronic programmable systems, including software, must be designed to ensure the repeatability, reliability and performance of these systems according to the intended use.	N/A			
6.2	Devices must be designed and manufactured in such a way as to minimise the risks of creating electromagnetic perturbation which could impair the operation of other devices or equipment in the usual environment.	N/A			

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
6.3	Devices must be designed and manufactured in such a way as to avoid, as far as possible, the risk of accidental electric shocks during normal use and in single fault condition, provided the devices are installed and maintained correctly.	N/A			
6.4 6.4.1	<p>Protection against mechanical and thermal risks</p> <p>Devices must be designed and manufactured in such a way as to protect the user against mechanical risks. Devices must be sufficiently stable under the foreseen operating conditions. They must be suitable to withstand stresses inherent in the foreseen working environment, and to retain this resistance during the expected life of the devices, subject to any inspection and maintenance requirements as indicated by the manufacturer.</p> <p>Where there are risks due to the presence of moving parts, risks due to break-up or detachment, or leakage of substances, then appropriate protection means must be incorporated.</p> <p>Any guards or other means included with the device to provide protection, in particular against moving parts, must be secure and must not interfere with access for the normal operation of the device, or restrict routine maintenance of the device as intended by the manufacturer.</p>	N/A			
6.4.2	Devices must be designed and manufactured in such a way as to reduce to the lowest possible level the risks arising from vibration generated by the devices, taking account of technical progress and of the means available for limiting vibrations, particularly at source, unless the vibrations are part of the specified performance.	N/A			
6.4.3	Devices must be designed and manufactured in such a way as to reduce as far as possible the risks arising from	N/A			

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	the noise emitted, taking account of technical progress and of the means available to reduce noise, particularly at source, unless the noise emitted is part of the specified performance.				
6.4.4	Terminals and connectors to electricity, gas or hydraulic and pneumatic energy supplies which the user has to handle must be designed and manufactured in such a way as to minimise all possible risks.	N/A			
6.4.5	Accessible parts of the devices (excluding the parts of areas intended to supply heat or reach given temperatures) and their surroundings must not attain potentially dangerous temperatures under normal use.	N/A			
7	Requirements for devices for self-testing Devices for self-testing must be designed and manufactured in such a way that they perform appropriately for their intended purpose taking into account the skills and the means available to users and the influence resulting from variation that can reasonably be anticipated in users' technique and environment. The information and instructions provided by the manufacturer should be easily understood and applied by the user.	N/A			
7.1	Devices for self-testing must be designed and manufactured in such a way as to: — ensure that the device is easy to use by the intended lay user at all stages of the procedure, and — reduce as far as practicable the risk of user error in the handling of the device and in the interpretation of the results.	N/A			
7.2	Devices for self-testing must, where reasonably possible, include user control, i.e. a procedure by which the user can verify that, at the time of use, the product will perform as intended.	N/A			

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
8 8.1	<p>Information supplied by the manufacturer</p> <p>Each device must be accompanied by the information needed to use it safely and properly, taking account of the training and knowledge of the potential users, and to identify the manufacturer.</p> <p>This information comprises the data on the label and in the instructions for use.</p> <p>As far as practicable and appropriate, the information needed to use the device safely and properly must be set out on the device itself and/or, where appropriate, on the sales packaging. If individual full labelling of each unit is not practicable, the information must be set out on the packaging and/or in the instructions for use supplied with one or more devices.</p> <p>Instructions for use must accompany or be included in the packaging of one or more devices.</p> <p>In duly justified and exceptional cases no such instructions for use are needed for a device if it can be used properly and safely without them.</p>	<p>Yes</p> <p>Yes</p> <p>Yes</p> <p>N/A</p>	<p>EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.</p> <p>EN ISO 18113-4:2011 In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling). Part 4: In vitro diagnostic reagents for self-testing</p>	<p>Product Dossier Section 5.2, Instructions for Use (IFU).</p> <p>Product Dossier Section 5.1, Labels.</p> <p>The "SARS-CoV-2 Antigen Packaging Operation Protocol" specified one IFU per test kit.</p>	<p>Management Protocol for Labels and IFU QM-00-001</p> <p>SARS-CoV-2 Antigen Packaging Operation Protocol PO1N4006</p>
8.2	<p>Where appropriate, the information to be supplied should take the form of symbols. Any symbol and identification colour used must conform to the harmonised standards. In areas for which no standards exist, the symbols and colour used must be described in the documentation supplied with the device.</p>	<p>Yes</p>	<p>EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.</p>	<p>Product Dossier Section 5.1, Labels.</p>	<p>Management Protocol for Labels and IFU QM-00-001</p>
8.3	<p>In the case of devices containing or a preparation which may be considered as being dangerous, taking account of the nature and quantity of its constituents and the form under which they are present, relevant danger symbols and labelling requirements of Directive 67/548/EEC and Directive 88/379/EEC shall apply. Where there is</p>	<p>N/A</p> <p>N/A</p>			

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	insufficient space to put all the information on the device itself or on its label, the relevant danger symbols shall be put on the label and the other information required by those Directives shall be given in the instructions for use. The provisions of the aforementioned Directives on the safety data sheet shall apply, unless all relevant information as appropriate is already made available by the instructions for use.				
8.4	The label must bear the following particulars which may take the form of symbols as appropriate: a) the name or trade name and address of the manufacturer. For devices imported into the Community with a view to their distribution in the Community, the label, the outer packaging, or the instructions for use shall contain in addition the name and address of the authorised representative of the manufacturer; b) the details strictly necessary for the user to uniquely identify the device and the contents of the packaging; c) where appropriate, the word 'STERILE' or a statement indicating any special microbiological state or state of cleanliness; d) the batch code, preceded by the word 'LOT', or the serial number; e) if necessary, an indication of the date by which the device or part of it should be used, in safety, without degradation of performance, expressed as the year, the month and, where relevant, the day, in that order; f) in case of devices for performance evaluation, the words 'for performance evaluation only'.	a) Yes b) Yes c) N/A d) N/A e) Yes f) Yes	EN ISO 18113-4:2011 In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling). Part 4: In vitro diagnostic reagents for self-testing EN ISO 15223-1:2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.	Product Dossier Section 5.2, Instructions for Use (IFU). Product Dossier Section 5.1, Labels. SARS-CoV-2 Antigen Packaging Operation Protocol Product Dossier Section 5.1, Labels (lot number specified). Product Dossier Section 5.1, Labels (shelf-life specified). Management Protocol for Labels and IFU	Management Protocol for Labels and IFU QM-00-001 SARS-CoV-2 Antigen Packaging Operation Protocol PO1N4006 BS-QA-007 Risk Management Control Procedures BS-QA-008 Design and Development Control Procedures

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	g) where appropriate, a statement indicating the <i>in vitro</i> use of the device; h) any particular storage and/or handling conditions; i) where applicable, any particular operating instructions; ¹ j) appropriate warnings and/or precautions to take; k) if the device is intended for self-testing, that fact must be clearly stated.	g) Yes h) Yes i) Yes j) Yes k) Yes		Intended use specified in IFU. Storage conditions specified in IFU and labels. Test procedures and relevant precautions specified in IFU. Relevant precautions and protective measures specified in IFU based on risk analysis. IFU specifies that the test is for self-testing.	
8.5	If the intended purpose of the device is not obvious to the user, the manufacturer must clearly state the intended purpose in the instructions for use and, if appropriate, on the label.	Yes	EN ISO 18113-4:2011 In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling). Part 4: In vitro diagnostic reagents for self-testing	Product Dossier Section 5.2, Instructions for Use (IFU). Product Dossier Section 5.1, Labels.	Management Protocol for Labels and IFU QM-00-001 BS-QA-008 Design and Development Control Procedures
8.6	Wherever reasonable and practicable, the devices and separate components must be identified, where appropriate in terms of batches, to allow all appropriate action to detect any potential risk posed by the devices and detachable components.	Yes	EN ISO 18113-4:2011 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 4: In vitro diagnostic reagents for self-testing	Product Dossier Section 5.1, Labels.	Lot Number Management Protocol PM-00-002

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
8.7	Where appropriate, the instructions for use must contain the following particulars:				
	a) the details referred to in section 8.4 with the exception of points d) and e);	a) Yes	EN ISO 18113-4:2011 In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling). Part 4: In vitro diagnostic reagents for self-testing	Product Dossier Section 5.2, Instructions for Use (IFU). Above contents included in IFU. Major raw materials and composition specified in IFU	Management Protocol for Labels and IFU QM-00-001
	b) composition of the reagent product by nature and amount or concentration of the active ingredient(s) of the reagent(s) or kit as well as a statement, where appropriate, that the device contains other ingredients which might influence the measurement;	b) Yes			BS-QA-008 Design and Development Control Procedures
	c) the storage conditions and shelf life following the first opening of the primary container, together with the storage conditions and stability of working reagents;	c) Yes	EN ISO 15223-1: 2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.	IFU specifies that the test card should be used immediately after opening.	
	d) the performances referred to in section 3 of part A;	d) Yes		Relevant information specified in IFU.	
	e) an indication of any special equipment required including information necessary for the identification of that special equipment for proper use;	e) N/A			
	f) the type of specimen to be used, any special conditions of collection, pre-treatment and, if necessary, storage conditions and instructions for the preparation of the patient;	f) Yes		Sample type, collection, and handling specified in IFU.	
	g) a detailed description of the procedure to be followed in using the device;	g) Yes		Detailed test procedures specified in IFU.	
	h) the measurement procedure to be followed with the device including as appropriate: — the principle of the method, — the specific analytical performance characteristics (e.g. sensitivity, specificity, accuracy, repeatability,	h) Yes		Test principle, product performance characteristics and sample handling methods specified in IFU. Testing does not require specific training if carried in accordance with the IFU.	

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	<p>reproducibility, limits of detection and measurement range, including information needed for the control of known relevant interferences), limitations of the method and information about the use of available reference measurement procedures and materials by the user,</p> <p>— the details of any further procedure or handling needed before the device can be used (for example, reconstitution, incubation, dilution, instrument checks, etc.),</p> <p>— the indication whether any particular training is required;</p> <p>i) the mathematical approach upon which the calculation of the analytical result is made;</p> <p>j) measures to be taken in the event of changes in the analytical performance of the device;</p> <p>k) information appropriate to users on:</p> <p>— internal quality control including specific validation procedures,</p> <p>— the traceability of the calibration of the device;</p> <p>l) the reference intervals for the quantities being determined, including a description of the appropriate reference population;</p> <p>m) if the device must be used in combination with or installed with or connected to other medical devices or equipment in order to operate as required for its intended purpose, sufficient details of its characteristics to identify the correct devices or equipment to use in</p>	<p>i) N/A</p> <p>j) N/A</p> <p>k) N/A</p> <p>—N/A</p> <p>—N/A</p> <p>l) N/A</p> <p>m) N/A</p>		<p>Product Dossier Section 5.1, Labels.</p>	

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	order to obtain a safe and proper combination;				
	n) all the information needed to verify whether the device is properly installed and can operate correctly and safely, plus details of the nature and frequency of the maintenance and calibration needed to ensure that the device operates properly and safely; information about safe waste disposal;	n) Yes		The reagent is the test, and does not require installing, maintenance and calibration.	
	o) details of any further treatment or handling needed before the device can be used (for example, sterilisation, final assembly, etc.);	o) N/A		Waste disposal specified in IFU.	
	p) the necessary instructions in the event of damage to the protective packaging and details of appropriate methods of re-sterilisation or decontamination;	p) N/A			
	q) if the device is reusable, information on the appropriate processes to allow reuse, including cleaning, disinfection, packaging and re-sterilisation or decontamination, and any restriction on the number of reuses;	q) N/A			
	r) precautions to be taken as regards exposure, in reasonably foreseeable environmental conditions, to magnetic fields, external electrical influences, electrostatic discharge, pressure or variations in pressure, acceleration, thermal ignition sources, etc.;	r) N/A			
	s) precautions to be taken against any special, unusual risks related to the use or disposal of the device including special protective measures; where the device includes substances of human or animal origin, attention must be drawn to their potential infectious nature;	s) N/A			
	t) specifications for devices for self-testing:	t) Yes		Above contents included in IFU.	

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 26 E-mail: 5 1 26
www.bosonbio.com

No.	Essential Requirement	Applicable ? (yes/no/na)	Method of Conformity	Evidence of Compliance and Location; Other remarks	SOP
	<p>— the results need to be expressed and presented in a way that is readily understood by a lay person; information needs to be provided with advice to the user on action to be taken (in case of positive, negative or indeterminate result) and on the possibility of false positive or false negative result, — specific particulars may be omitted provided that the other information supplied by the manufacturer is sufficient to enable the user to use the device and to understand the result(s) produced by the device, — the information provided must include a statement clearly directing that the user should not take any decision of medical relevance without first consulting his or her medical practitioner, — the information must also specify that when the device for self-testing is used for the monitoring of an existing disease, the patient should only adapt the treatment if he has received the appropriate training to do so; date of issue or latest revision of the instructions for use.</p>				

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China
Tel : 5 1 2e E-mail: 5 1 2e
www.bosonbio.com

9 RISK MANAGEMENT

Risk Analysis and Control Summary

In accordance with the standard requirements of EN ISO 14971:2012, the R&D department has organized relevant personnel to conduct risk analysis using the failure mode and effects analysis (FMEA) on the following aspects:

- 1) Risk analysis on the users, patients, manufacturing risks, etc. during the input phase of design development;
- 2) Risk analysis of the production processes (suppliers, material procurement, process control, process safety information) and other aspects during the trial production phase;
- 3) Risk analysis of the labels and IFU;
- 4) Risk analysis of transport, storage, and customer use before market launch.

Through risk analysis and assessment, a series of risk control measures were established and validated by the R&D department, including:

- 1) study on common interfering substances to SARS-CoV-2 virus to avoid potential effects on the test results;
- 2) study on the storage, transport and in-use stabilities within the product shelf-life to avoid effect of stability on the test results;
- 3) verification of the reagent processes to ensure that the production processes are reasonable and feasible.
- 4) Inspection of every batch of reagents, verification of maximum batch quantity, and specification of maximum batch quantity in BS-QA-011 "Production and Service Control Procedures" to avoid within-lot and lot-to-lot variance;
- 5) Tight sealing in product design and specification in the IFU on how to handle reagents after use to avoid the impact of testing on user safety and the environment;
- 6) Illustrations and accompanying texts in the IFU to specify the operation procedures and precautions to avoid testing errors in private use that could lead to erroneous test results.

Measures have been taken to minimize the risks that may arise in each process. See "Appendix 1 – Risk Management Report" for more comprehensive risk analysis and control measures. Based on analysis of the residual risks and the risk guidelines, the medical benefits of the intended use were evaluated to outweigh the residual risks.

**Xiamen Boson Biotech Co.,Ltd.**

90-94 Tianfeng Road, Jimei North Industrial Park,

Xiamen, Fujian 361021, P.R. China

Tel : [REDACTED] E-mail: [REDACTED]

www.bosonbio.com

Finally, the "Post-Marketing Surveillance Control Procedures" was established to specify the responsibilities of each department to collect information on production and post-production activities as required, and to provide continuous follow-up on the post-marketing status of the SARS-CoV-2 antigen test kits. See "Appendix 2 - Post-Marketing Surveillance Control Procedures" for the complete document.

Xiamen Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian 361021, P.R. China

Tel : 5.1.2e E-mail: 5.1.2e
www.bosonbio.com

10 LIST OF APPENDICES

In numerous places in this Technical File, references are made to the following appendices which are to be regarded as binding parts of this document.

Appendix 1 - Risk Management Report.pdf

Appendix 2 - Post-Market Surveillance Control Procedures

Appendix 3-1 - Analytical Sensitivity Studies.pdf

Appendix 3-2 - Analytical Sensitivity Studies - Nasal Sample Validation Report.pdf

Appendix 4 - Analytical Specificity - Cross Reactivity Studies.pdf

Appendix 5 - Analytical Specificity - Interference Studies.pdf

Appendix 6 - Precision Studies.pdf

Appendix 7 - Hook Effect Studies.pdf

Appendix 8 - Storage Stability Studies.pdf

Appendix 9 - Freeze-Thaw Stability Studies.pdf

Appendix 10 - Transport Stability Studies.pdf

Appendix 11 - Heat Stability Studies.pdf

Appendix 12 - Open Pack Stability Studies.pdf

Appendix 13 - Clinical Evaluation Report.pdf

Appendix 14 - Private-Use Usability Study.pdf