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Aanvrager: Tintele Global Co, Ltd - March 28, 2021

Handelsnaam: **VI-PASS**
SARS-CoV-2 antigeen rapid test kit for home use

Producent: BEIJING LEPU MEDICAL TECHNOLOGY CO, LTD
Building 7-1 no. 37 Chaoqian Road, Changing district, Beijing,
102200, PRChina

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I. Principe



- Card-like format, detects the nucleocapsid (N) protein of the virus.
- Gold-labelled SARS-CoV-2 N protein monoclonal antibodies immobilized on the test area, with corresponding antibodies in the quality control area.
- Monster/specimen: Anterior Swabbed specimens from the nasopharynx (2-3 cm inside nose)
- result in 15 minutes.

*subject to change according to EU legislation and manufacturer's design


2019-nCoV Antigen Rapid Test Kit (Colloidal Gold Immunochromatography)


Home Test

Operating Steps

- Wash your hand, take out test card from outer package.
 
- Place test card flat on table, remove cover-layer of adhesive.
 
- Take out swab from stick-end, refer to standard nasal swab specimen collection to collect specimen.
 

Note: Do not touch the swab head.

Note: Sampling in both nasal cavity sample is required.
- Insert the swab head into well A from the bottom of well B.
 

Keep the card flat on table.
- Add 6 drops of the 2019-nCoV antigen test sample treatment solution to well A, then rotate 2 rounds, each direction in the buffer.
 

Keep the card flat on table.

Note: False negative results may occur if the sample swab is not turned before closing the test card.

Note: Do not rotate the swab while dropping the sample.
- Fold the left side over, fit two sides together completely, start timing.
 

Keep the card flat on table.
- Wait for the appearance of purple-red line. Test results should be read within 15-20 minutes.
 

Keep the card flat on table, Do not move the test card.

15-20 min

Note: False results can occur if the card is disturbed/moved.

Note: False results can occur if the test results are read before 15 minutes or over 20 minutes.
- After test, put the test card, swab, and dilution bottle into outer package and seal it tightly. Dispose the bag in waste container according to local laws and regulations.
 


2019-nCoV Antigen Rapid Test Kit (Colloidal Gold Immunochromatography)

Home Test

Interpretation of Test Results

- Positive (+): A purple-red band appears in the Control Line (C) and Test Line (T).

Solid Line

Control Line
Test Line
Positive


Faint Line

Control Line
Test Line
Positive

OR


The bottom line can be very faint Any pink/purple line visible here is positive.
- Negative (-): Only the Control Line (C) shows a purple-red band. No purple-red band appears in the Test Line (T).

Control Line
Test Line
Negative
- Invalid: If "no purple-red band appears in Control Line (C)" and "a blue band appears in the Control Line (C)", it indicates that the operation process is incorrect or the test paper has been damaged. In this case, please read the instruction manual carefully again and retest with a new test paper. If the problem persists, please stop using this batch of products immediately and contact your local supplier.


Control Line
Test Line
Invalid

Control Line
Test Line
Invalid

Control Line
Test Line
Invalid

Control Line
Test Line
Invalid

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- Step 1:** Use a swab to collect anterior nasal specimen (2-3cm in nose)
- Step 2:** Peel off adhesive.
- Step 3:** Insert swab from well B to well A.
- Step 4:** Add 6 drops of buffer to well A rotate the shaft, two rounds each direction.
- Step 5:** Fold the left side over, and wait for 15 minutes to see the result appear on the front of the card.

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(appendix II)

		Checklist Essential Requirements IVDD		001/08/2020	
				ID: 5982	
				Customer Number 5982	
Company BEIJING LEPU MEDICAL TECHNOLOGY Co., Ltd					
Products/Product groups SARS-CoV-2 Antigen Rapid Test kit for self use - VI PASS					
Conformity assessment procedure section 2 to 5 in annex III of IVDD 98/79/EC			Categorization IVD others		
ER	Applicable standards, other documents	Documented evidence (test reports, records, others)	NA*	fulfilled	
I. General requirements					
1	The devices must be designed and manufactured in such a way that, when used under the conditions and for the purposes intended, they will not compromise, directly or indirectly, the clinical condition or the safety of the patients, the safety or health of users or, where applicable, other persons, or the safety of property. Any risks which may be associated with their use must be acceptable when weighed against the benefits to the patient and be compatible with a high level of protection of health and safety.		<input type="checkbox"/>	<input checked="" type="checkbox"/>	
2	The solutions adopted by the manufacturer for the design and construction of the devices must conform to safety principles, taking account of the generally acknowledged state of the art. In selecting the most appropriate solutions, the manufacturer must apply the following principles in the following order: – eliminate or reduce risks as far as possible (inherently safe design and construction); – where appropriate take adequate protection measures in relation to risks that cannot be eliminated; – inform users of the residual risks due to any shortcomings of the protection measures adopted.		<input type="checkbox"/>	<input checked="" type="checkbox"/>	
3	The devices must be designed and manufactured in such a way that they are suitable for the purposes referred to in Article 1(2)(b), as specified by the manufacturer, taking account of the generally acknowledged state of the art. They must achieve the performances, in particular, where appropriate, in terms of analytical sensitivity, diagnostic sensitivity, analytical specificity, diagnostic specificity, accuracy, repeatability, reproducibility, including control of known relevant interference, and limits of detection, stated by the manufacturer. The traceability of values assigned to calibrators and/or control materials must be assured through available reference measurement procedures and/or available reference materials of a higher order.		<input type="checkbox"/>	<input checked="" type="checkbox"/>	
4	The characteristics and performances referred to in sections 1 and 3 must not be adversely affected to such a degree that the health or the safety of the patient or the user and, where applicable, of other persons, are compromised during the lifetime of the device as indicated by the manufacturer, when the device is subjected to the stresses which can occur during normal conditions of use. When no lifetime is stated, the same applies for the lifetime reasonably to be expected of a device of that kind, having regard to the intended purpose and the anticipated use of the device.		<input type="checkbox"/>	<input checked="" type="checkbox"/>	
5	The devices must be designed, manufactured and packed in such a way that their characteristics and performances during their intended use will not be adversely affected under storage and transport conditions (temperature, humidity, etc.) taking account of the instructions and information provided by the manufacturer.		<input type="checkbox"/>	<input checked="" type="checkbox"/>	

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Clinical performance

A total of 508 clinical specimens based on nucleic acid assay (PCR) were collected, including 243 positive specimens and 265 negative specimens. After comparing this product with nucleic acid assay (PCR) through the collected clinical samples, the results are summarized as follows:

SARS-COV2 Antigen Rapid Test Kit	Nucleic acid assay (PCR)	
	Positive	Negative
Positive	231	1
Negative	12	264
Analysis of sensitivity	95.06% (95%CI:91.57%~97.15%)	/
Analysis of specificity	/	99.62% (95%CI:97.89%~99.93%)

Performance against the Comparator Method-by Cycle Threshold Counts.

SARS-COV2 Antigen Rapid Test Kit	Nucleic acid assay (PCR)	
	Positive(Ct≤32)	Positive(Ct≤25)
Positive	227	202
Negative	8	3
sensitivity	96.60% (95%CI:93.43%~98.27%)	98.54% (95%CI:95.79%~99.50%)

Foto 10

II. Validatiestudies

1. analysis Specificity Investigation Report (appendix III)

by Beijing Lepu Medical Technology Co., LTD

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Contents

1 Overview	2
2 Information about the kit used in this investigation.....	2
3 Investigation method	2
3.1 Content of cross-reaction validation.....	3
3.2 Validation information.....	3
3.3 Validation information.....	3
4 Investigational results.....	4
4.1 Test results of nasal swab samples from 20 healthy subjects.....	4
4.2 Test results of pathogen samples and recombinant N protein.....	5
4.3 Test results of interferent.....	8
5 Conclusion.....	12

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1. **Stability study report (appendix VI)**
by Beijing Lepu Medical Technology Co., LTD

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- Directive: 98/79/EC of the European Parliament
Section 2 to 5 in Annex III of IVDD 98/79/EC

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V. EU common list registratie SARS-CoV-2 Antigeen Rapid test for professional use (appendix XII)

- link: https://ec.europa.eu/health/sites/health/files/preparedness_response/docs/covid_19_rat_common_list_en.pdf

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