

<b>HQTS</b>	<b>INSPECTION REPORT</b> Pre-Shipment Inspection	Report No: H-20062555
-------------	---	-----------------------

INSPECTION INFORMATION			
Client	MEDIQ		
Supplier	Hunan Sanrui Biotchnology Co., Ltd		
Factory	Hunan Sanrui Biotchnology Co., Ltd		
Inspection Instruction	GWI-H-004-HQT		
P.O. No.	5.1.1c	P.O. Quantity	3600000 Pieces
Item No.	SRNS-SM-A	Sampling Plan	ISO2859 : Single, Normal
Inspection Date	09/06/2020 to 10/06/2020	Inspection Level	G-II
Inspection Sequence	Original inspection	Sample Size	800 Pieces
Inspection Location	ChangSha,China	AQL	Critical:0,Major:1.5,Minor:4.0
Product Description	SURGICAL MASK TYPE IIR		



INSPECTION RESULT SUMMARY			
Category	Conform	Not Conform	Pending
1.Quantity	X		
2.Workmanship	X		
3.Style, Color & Documentation	X		
4.Field Test & Measurement	X		
5.Shipping Mark & Packaging			X
<b>OVERALL CONCLUSION</b>	<b>PENDING</b> for client's evaluation		

SPECIAL ATTENTION POINT (SAP)	
Approved by HQTS Supervisor:	5.1.2e

*This report reflects the facts as recorded by HQTS at the time and place of inspection. It does not relieve the manufacturers from their contractual obligations nor prejudice client's right for compensation for any apparent and/or hidden defects not detected during our random inspection or occurring thereafter.*

Report No: H-20062555

1	Quantity	NIL
2	Workmanship	NIL
3	Style, Color & Documentation	NIL
4	Field Test & Measurement	NIL
5	Shipping Mark & Packaging	
5.1	The lot No of label was not stuck on the export carton during the inspection, but the factory claimed it would be done before shipment.	







#### INSPECTION DETAILS

1. Quantity					Conform			
P.O. No.	Item No.	P.O. Qty	Shipment Qty		Packed Qty		Sample Size	
		(Piece)	(Piece)	(Carton)	(Piece)	(Carton)	(Piece)	(Carton)
5.1.1c	SRNS-SM-A	3600000	3600000	1440	3490000	1396	800	38
Total		3600000	3600000	1440	3490000	1396	800	38

2. Workmanship		Conform		
Defects Description		Critical	Major	Minor
Item#: SRNS-SM-A				
Damaged on ear loop			2	
Obvious different length of ear tape(1-2cm)			3	
Slightly different length of ear tape(0.5-1cm)				1
Uneven rolling edge at surface(1cm)				3
Total Found:		0	5	4
Allowed:		0	21	21
Sample Size:		800	800	315

The following field test included in workmanship check

Defect Photos

	
Defect-damaged on ear loop-ma	Defect-damaged on ear loop-ma
	
Defect-Obvious different length of ear tape (1-2cm)-ma	Defect-Obvious different length of ear tape (1-2cm)-ma
	
Defect-Obvious different length of ear tape(1-2cm)-ma	Defect- Slightly different length of ear tape(0.5-1cm)

<p>Defect-uneven rolling edge at surface(1cm)-mi</p>	<p>Defect-uneven rolling edge at surface(1cm)-mi</p>
	<p>NA</p>
<p>Defect-Uneven rolling edge at surface(1cm) -mi</p>	

3. Style, Color & Documentation	<b>Conform</b>		
DESCRIPTION	Conform	Not Conform	Pending
The inspected samples were consistent with	-	-	-
Purchase Order	X		
Digital photo	X		
Instruction manual and Document Verification	-	-	-
Client Comment	-	-	-

4. Field Test & Measurement			<b>Pending</b>		
Item	Specification	Sample Size	Conform	Not Conform	Pending
4.1	Wear check	5 samples	X		
4.2	Product dimension check	5 samples	X		
4.3	Product net weight check	5 samples	X		
4.4	Smell check	5 samples	X		
4.5	Water leakage test	5 samples	X		
4.6	Rub test	5 samples	X		
4.7	Quantity per box check	All selected	X		

Report No: H-20062555

			cartons			
4.8	Carton dimension check		1 carton	X		
4.9	Carton gross weight check		1 carton	X		
4.10	Carton drop check		1 carton	X		
4.11	Tension test on ear string		5 samples	X		
4.12	Bacterial filtration efficiency test		5 samples	X		
4.13	Synthetic blood penetration test		5 samples	X		
4.14	Differential pressure test		5 samples	X		
4.14	Microbial cleanliness test		5 samples	X		
4.14	Fire proof to melt-blown fabric test		5 samples	X		

Item No.	Carton dimension (cm) ( $\pm 5\%$ )		Carton gross weight (kg)/ ( $\pm 10\%$ )	
	Specification	Result	Specification	Result
SRNS-SM-A	57. 5x40x49	60.5x41.3x50.2	11	11.54

Item No.	Unit package dimension(cm)		Unit package weight (g)	
	Specification	Result	Specification	Result
SRNS-SM-A	N/A	19.0x11.0x9.0, 19.0x11.0x9.0, 19.0x11.0x9.0, 19.0x11.0x9.0, 19.0x11.0x9.0	N/A	220, 225, 221, 225, 222

Item No.	Product dimension(cm)		Product net weight (g)	
	Specification	Result	Specification	Result
SRNS-SM-A	18x10	17.5x9.5, 17.5x9.5, 17.5x9.5, 17.5x9.5, 17.5x9.5	N/A	2.9, 2.8, 2.9, 2.8 ,2.9

Item No.	Fabric weight (gsm)	
	Specification	Result
SRNS-SM-A	-	The first layer of cloth:25 gsm The middle layer of cloth:27 gsm The third layer of cloth:29 gsm

<b>5. Shipping Mark &amp; Packaging</b>	<b>Pending</b>
5.1 Shipping mark photos	

Report No: H-20062555



## 5.2 Packaging photos



## Packaging Description:

- 50 pieces of products were packed into a polybag Choose.
- 1 piece of polybag was packed into an color box.
- 50 color boxes were packed into a 5-ply corrugated export carton.
- Export carton was sealed with gummed tape.


## 6. Reference Note

6.1	For microbial cleanliness test and bacterial filtration efficiency test was not conducted, since the test need 48 hours of bacterial and microbial culture during the inspection, Our inspector only pull out the products, The test process and results were inspected and would be provided later by the factory.
6.2	All factory laboratory experiments, because they are sophisticated experiments, Our inspector only witnessed the experiment and result was just for reference.
6.3	Carton number was not assigned by the factory yet.
6.4	Manufacturer did not allow us to draw any production sample for reference.
6.5	Our inspector required the factory to repair all defective samples before shipment.
6.6	Upon arrival, all products were finished production, among which 3490000 pieces (1440 cartons) were packed, and the balance was under packing.


7. Photos




Warehouse




Test for strength of belt




Test for strength of belt




Test for strength of belt report



Synthetic blood penetration test



Synthetic blood penetration test



Synthetic blood penetration test

Synthetic blood penetration test

Differential pressure test

Differential pressure test

Differential pressure test

Fire proof to melt-blown fabric test


欢迎使用  
广州际包装设备有限公司  
GBN701口罩透气阻力及压力差测试报告

执行标准:  
实验环境:  
材料编号:  
实 验 员:  
实验类型:  
实验单位:  
实验日期: 2020-06-09

设定流量(L/min): 8.00  
设定时间(s): 10  
测试面积(cm<sup>2</sup>): 4.00


压方差 Pa/(cm<sup>2</sup>)

实验(1): 33.1  
实验(2): 26.7  
实验(3): 23.6  
实验(4): 21.0  
实验(5): 20.8  
平均值: 25.0



Fire proof to melt-blown fabric test


Material structure check



Particulate matter filtration efficiency test



Particulate matter filtration efficiency test



Particulate matter filtration efficiency test report




Wear check

	
<p>Export carton</p>	<p>Main mark</p>
	
<p>Main mark</p>	<p>Side mark</p>
	
<p>Main mark</p>	<p>Side mark</p>




Packing view


Packing view




Specification




Specification



Specification



Specification



Specification

Specification




Color box



Color box








Color box



Color box

 <p>The image shows the exterior of a white and blue cardboard box for 'Medical surgical masks (non-sterile)'. It features a blue surgical mask illustration, a 'G.C. PASSED' logo, and a '50' count indicator.</p>	 <p>The image shows the back of the box with detailed text in Chinese and English, including 'Color box', 'Contents list', 'Instructions', and 'Precautions'.</p>
<p>Color box</p>	<p>Color box</p>
 <p>The image shows the box from an internal perspective, revealing the white protective paper lining and the organized arrangement of the masks inside.</p>	 <p>The image shows the box open, with a clear plastic bag containing a roll of blue surgical masks and a small instruction card placed next to it.</p>
<p>Packing view</p>	<p>Packing view</p>
 <p>The image shows a 'Certificate of qualification' (产品合格证) for the masks. It includes fields for product name, specifications, standards, and manufacturer information, with a red circular stamp.</p>	 <p>The image shows a single blue surgical mask in its clear plastic packaging, highlighting the three layers and the ear loops.</p>
<p>Certificate of qualification</p>	<p>Packing view</p>

	
<p>Product view</p>	<p>Product view</p>
	
<p>Product view</p>	<p>Product view</p>
	<p>NA</p>
<p>Product view</p>	

Report No: H-20062555

<b>Inspector:</b>	2250	<b>Auditor:</b>	666
		<b>End of Report</b>	_____