



- Fast
- Cost-Effective
- Easy to Use
- Easy to transport
- Accurate



Edinburgh Genetics

ActivXpress+ COVID-19 Antigen Complete Test Kit

Edinburgh Genetics ActivXpress+ COVID-19 Antigen Complete Testing Kit enables reliable Sars-Cov-2 antigen testing outside laboratories in 15 minutes.



Benefits

Fast

Results are available in 15 minutes, allowing rapid isolation of patients, minimising potential infections.

Cost-Effective

Testing can be done without expensive laboratory equipment.





Benefits


Easy to Use

Staff at the decentralised locations including primary health care and community-level health care can administer the test with all equipment required provided. Results can be interpreted visually without the need for an expensive reader.





Benefits



Easy to transport

The product can be transported and stored at room temperature, enabling reliable and cost-effective testing in the world's most inaccessible places.



Benefits

Accurate

Our reliable tests produce market-leading results on sensitivity and specificity that exceed the standards set by the WHO, without cross-reactivity to seasonal coronaviruses.





Packaging

The Edinburgh Genetics ActivXpress+ COVID-19 Antigen Complete Testing Kit comes as a box of 1, 10 or 20



Individual Pack Box



Multipack Box



Packaging

Catalogue Numbers

EGCV0101

1 x Test Cassette, 1 x Sterilised nasopharyngeal swab, 1 x reagent in tube with dropper

EGCV0101A

10 x Test Cassette, 10 x Sterilised nasopharyngeal swab, 10 x reagent in tube with dropper

EGCV0101B

20 x Test Cassette, 20 x Sterilised nasopharyngeal swab, 20 x reagent in tube with dropper

EGCV0101M


1 x Test Cassette, 1 x Sterilised oropharyngeal swab, 1 x reagent in tube with dropper

EGCV0101MA

10 x Test Cassette, 10 x Sterilised oropharyngeal swab, 10 x reagent in tube with dropper

EGCV0101MB

20 x Test Cassette, 20 x Sterilised oropharyngeal swab, 20 x reagent in tube with dropper





About Edinburgh Genetics

As a ISO 13485:2006 certified manufacturer, we have been developing high quality, innovative and accessible IVD devices that are used clinically in countries around the world.





Contact Us



Contact Us

www.eggenetics.com

info@eggenetics.com

(44) 131 261 6686

Glencorse Building, Pentlands Science Park,
Penicuik, Midlothian, UK, EH26 0PZ